



MODEL NO.: YSI-155

YSI ELISA PLATE ROTATOR

For Elisa plate shaking at adjustable speed from minimum to maximum required for Elisa test. Provided with digital speed controller cum indicator. The Rotator having platform of size 300 x 300 which can take Elisa plates in the special framed cavities. An imported timer 0-1 hr. is also fitted to control the shaking duration. Supplied complete with cord & plug. Suitable to work on 220 V Single Phase 50 HZ Ac supply.



MODEL NO.: YSI-156

YSI PLASMA SEPARATOR

For efficient, Clear Separation of Blood Components. The stand is designed to Extract Blood Components from a Blood Collecting Bag to a Satellite Bag or from a Satellite bag to a transfer bag. Keeps constant gentle pressure for separating plasma with safe transfer. Blood bag is kept in between pressure plate and the uniform wall of the equipment. After releasing the pressure it is easy to transfer from blood bag to satellite bag. It is easy to operate. Accommodates bags upto Capacity : 500 ml



MODEL NO.: YSI-158

YSI SYRINGE NEEDLE DESTROYER

Scientific research has discovered that used syringes and needles are very dangerous to humans that can cause incurable diseases like AIDS, Hepatitis and many more. A unique equipment has been developed by YSI to destroy both i.e. needle and syringes in a Safe, Eco-friendly and Cost effective manner.

After switching ON, insert the used needle to the small hole provided and within seconds the needle will be destroyed. Insert the syringe on the top of the machine in the big 2ml hole provided with a cutting blade, just pull the lever and syringe will be cut into two. A tray is built in to collect Needle residue & Syringe stubs. Specially designed electrodes generate high temperature to instantly incinerate the needle.



MODEL NO.: YSI-158A

YSI SYRINGE NEEDLE DESTROYER

Same as ..YSI-158 but provided with Battery back up With built in power source trickle charger. Battery is constantly charged when main in ON. It shifts automatically to battery operation in case of power failure, without disturbance in the operation An overload trip device is also provided. Beside the needle residue tray, a syringe container (Optional) is provided to collect syringe bodies after cutting.